

Vertex Dental

Vertex™ Quint Teeth Mould Chart



Vertex™ Quint Synthetic Resin Teeth Upper Anteriors - Square




Vertex™ Quint Synthetic Resin Teeth Upper Anteriors - Ovoid




Advantages of the Vertex™ Quint Teeth:

- Excellent bonding characteristics with Vertex™ Denture Base Acrylics
- Wide range of models
- Available in A1 until D4 colour classification
- Innovative 5 layer concept
- Easy to grind with cervical consistent color ceramic effect
- Double cross linked and a high molecular weight






Vertex™ Quint Synthetic Resin Teeth Upper Anteriors - Tapering





Vertex™ Quint Synthetic Resin Teeth Lower Anteriors





Vertex™ Quint Synthetic Resin Teeth Posteriors





Posterior measurements



All dimensions are in millimeters

Vertex™ Quint Synthetic Resin Teeth - Combined Articulations

Anteriors		Posteriors			
Upper	Lower	A	B	C	D
M1	AL1	8	12	11	14
M2	AL2	8	13	11	14
M3	AL4	7	13		15
M4	AL11	9	7		15
M5	AL8/AL9	10	9		15
R1	AL1	8	12	11	14
R2	AL2	7	13	11	14
R3	AL5/AL5	9	13	11	15
R4	AL6	9	7		15
R5	AL9/AL10	10	9		15
R6	AL7/AL11	9	7		15
R7	AL3	8	13	11	14
R8	AL11	7	9	13	15
R9	AL12/AL11	7	8	13	14
J1	AL1	8	11	12	14
J2	AL2	8	11	12	14
J3	AL12	7	9	11	14
J4	AL11/AL12	9	7	10	15
J5	AL9/AL10	10	9	7	15
J6	AL7/AL11	10	9	7	15

Vertex™ Quint Synthetic Resin Teeth - Shades

Available in 16 shades and tones:

A1, A2, A3, A3.5, A4

B1, B2, B3, B4

C1, C2, C3, C4

D2, D3, D4

The dimensions are in millimeters.

Vertex™ Quint Teeth meet all the requirements of EN-ISO 22112:2006 and are not fluorescent.

Vertex Dental

Vertex-Dental B.V.

Johan van Oldenbarneveltlaan 62 | 3705 HJ Zeist
 P.O. Box 10 | 3700 AA Zeist | The Netherlands
 Tel. +31 (0)30-69 76 749 | Fax +31 (0)30-69 55 188
 info@vertex-dental.com | www.vertex-dental.com

BAT1306

